

4/29/17

TITLE

PROJECT

Continued from page


Inoculating P_{int} constructs in LB

	P _{int} -WT	P _{int} -CTD	
O/N	4.685	5.865	(corrected w/dilution)
Add to 1L	10.7 mL	8.5 mL	
8:54a	0.050	0.045	
9:55a	0.138	0.155	
10:31a	0.364	0.366	
10:46a	0.503	0.507	← Induced 1mM IPTG
5:22	4.86	2.74	← harvest

Spun down cells, I and U sample. Each had OD ≈ 0.5 and 1 mL cells will resuspend in Urea (6M) / 20 mM Tris / 50 mM NaCl. Some volume left over in LB, I'll add 70 µL, then dye.

Making my own stock of loading dye. Lab stock is 4x non-reducing. I need to add 15 µL BME to 1 mL dye.

2.86 g of CTD cells, 3.43 g of WT cells.



Note: Samples were mixed 1:1 w/ 2x LDS dye + BME and heated to 90 °C for 5 min before loading. Not viscous @ all. Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY


DATE

PROPRIETARY INFORMATION

TITLE PROJECT

Continued from page

5 Gel lanes are on prev. page. Pattern is CTD
 U-I, WT U-I CTD
 10 is measuring 14 kDa,
 FL WT pin1 is around
 15 ~~10~~ kDa. This is pretty
 close to the expected
 values (see next page.)



5/1/2017

Lysis Buffer / Wash Buffer (1L)			Stock
150 mM NaCl → 30 mL	30 mL		5M
20 mM PIPES ^{NEPES} (pH 7.5) →	20 mL		1M
20 mM Imidazole (pH 7.5) →	5.7 mL		3.5M
2 mM DTT	0.31 g (solid)		10M

Wash Buffer Elution Buffer (500mL)		
150 mM NaCl → ¹⁵ 30 mL		
20 mM NEPES → 10 mL		
800 mM Imid → 115 mL		
2 mM DTT → 0.16 g		

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

Continued from page

Add ~40 mL lysis buffer to each pellet (+ 0.5 mg/mL lysozyme). Vortex ~~for 10 min~~ and shake for 10' to resuspend pellet in each conical tube. Now rocking on ice for another 10'. Sonicate 6' total, 2' pulse, 1' rest at power 6. Kept on ice during sonication, though felt warm. Spin at 19,000 rpm in JA 25.50 for 45'. Took sample for gel.


Column turned black w/ 2 mM DTT. Running 4 mL min⁻¹, reus HisTrapp from Dinusha. Starting w/ CTD.

$$60 \text{ mL} \frac{\text{min}}{4 \text{ mL}} = 15 \text{ min gradient from } 0-100\% \text{ B.}$$

Collecting 5 mL fractions

Result for CTD is below I will pool A3-A6

UNICORN 5.20 (Build 500)
Result file: c:\...\dinusha\20170501 Unlabeled Pin1 CTD HisTrap



TITLE

PROJECT


DATE

Continued from page

Running WT, same settings, different column. Shape of elution is odd. Bad column? Overloaded? I will pool A6-B12.

UNICORN 5.20 (Build 500)

Result file: c:\...\dinusha\20170501 Unlabeled Pin1 WT Histrap



Need dialysis buffer. Will use 2L of

150 mM NaCl		60 mL
20 mM HEPES	→ pH 7.5	40 mL
2 mM DTT		0.62 g

WT Pin1 pooled sample (1mm=1)	$A_{280} = 0.373 = 177 \mu\text{M}$
CTD " " " "	$A_{280} = 0.277 = 396 \mu\text{M}$

$$\frac{177 \times 10^{-6} \text{ mol}}{\text{L}} \cdot \frac{20963 \text{ g}}{\text{mol}} \times \frac{1000 \text{ mg}}{\text{g}} \times 0.04 \text{ L} = 149 \text{ mg WT Pin1}$$

$$\frac{396 \times 10^{-6} \text{ mol}}{\text{L}} \cdot \frac{15039 \text{ g}}{\text{mol}} \cdot \frac{1000 \text{ mg}}{\text{g}} \cdot 0.02 \text{ L} = 120 \text{ mg CTD Pin1}$$

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

DATE

Continued from page

Prep recommends 10/mg. Stock is 400/ μ L. 3.6 μ L for
5 WT Pinl and 3 μ L for CTD. Stir Stirred O/N @ 4degC
in Emerson fridge (4:30 pm)

10 Ran 2 gels, one for each prep

15 WT
M/Induced/Lysed/Supernatant/FT/Wash/Pool/Blank/M

20 CTD
M/Induced/Lysed/Supernatant/FT/Wash/Pool/M

25 Gels are asymmetric, so I should be able to tell them
apart easily.

30 5/2/17

Removed from dialysis (looked good) and put in fridge @
8:30 am.

35 ~~Add back to~~ Took 2x1ml of benz. beads
and added dialysis buffer to equilibrate. Rocked
40 samples w/ beads ~5 min. Spinning down. 4k ref
for 5 min. Will filter in 0.45 μ m filter.

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT


5/2/17 DATE

Continued from page

Eluting CTD through column. Running 5 mL/min, after load switch to 2% B to help w/ nonspecific binding. That led to a 2nd peak, which I'll keep before pooling. Will pool A2-A9, run A8 separately (but still pool it)


UNICORN 5.20 (Build 500)

Result file: c:\...\dinusha\20170502 Unlabeled Pin1 CTD His Tag Removal



No more Ni-NTA - I will add 10 mM DTT to pooled protein to protect it during conc, gel filt. (400 µL)

Some thing w/ WT. I will set aside some ~~A11~~ ~~A12~~ ^{B11} fraction & run separately. Will pool A2-B12. More volume here, will still add 10 mM DTT (600 µL)



Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE


PROJECT

DATE

Continued from page

UNICORN 5.20 (Build 500)

Result file: c:\...\dinusha\20170502 Unlabeled Pin1 WT His Tag Removal



Concentrations

	Volume	C (µM)	mg/ml	µmol/mol	mg
CTD 1mm	40 mL	147	1.93	5.58	77
WT Pin1 #1 (A2-A7) 1mm	30 mL	70	1.30	2.1	40
WT Pin1 #2 (A8-B11) 1mm	35 mL	108	2.00	3.78	70

$$\frac{x \text{ } \mu\text{mol}}{L} \cdot \frac{1 \text{ mol}}{10^6 \text{ } \mu\text{mol}} \cdot \frac{\text{Mg}}{\text{mol}} \cdot \frac{1 \text{ mg}}{0.001 \text{ g}} \cdot \frac{1 \text{ L}}{1000 \text{ mL}} = \frac{x \cdot M}{10^6} \frac{\text{mg}}{\text{mL}}$$

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

DATE

Continued from page

Need to make gel filtration buffer. Given application (activity) I will use the following

Total Vol: 1.5 L



	<u>Stock</u>	<u>Vol</u>
50 mM NaCl	5M	15 mL
10 mM Imidazole pH 6.5	3.5M	4.3 mL
1 mM TCEP	0.5M	3 mL
0.02% (w/v) NaI₂	—	0.3 g

Forgot will add later.

Gels from yesterday are below. See p. 46 for lanes.

WT Pnl

CTD Pnl



Note strong 15 kDa band in both WT and CTD growths while the FT looks like CTD lost protein, the presence in WT means that it's not likely.

SIGNATURE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

DATE

Continued from page


Ran concentration of CTD sample - 0.632 for 1mm. That's
Volume is approx 5 mL. Conc is 904 μ M.

Conc of WT is approx since UV is so high. Nanodrop
registers 1.997 @ 1mm, but I think that's probably
not reliable. As a working value, the cis approx 952 μ M
@ A 1:5 dilution (2:10 μ l) gives A=0.453, so its close.

UNICORN 5.20 (Build 500)

Result file: d:\...\dinusha\20170502 Pin1 CTD Superdex 75 26 60001

- 20170502 Pin1 CTD Superdex 75 26 60001:10 UV1 280nm
- 20170502 Pin1 CTD Superdex 75 26 60001:10 UV2 260nm
- 20170502 Pin1 CTD Superdex 75 26 60001:10 Cond
- 20170502 Pin1 CTD Superdex 75 26 60001:10 Fractions



I will pool C2-C8. Adding NaN₃ to 0.02%.

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

5/3/2017 DATE

Continued from page

$$10\% \frac{(35 \times 0.02)}{10\%} = 70 \mu\text{L of } 10\% \text{ NaN}_3$$

Azss of CTD is 0.093 (1mm). Conc is 133 μM ($V=35\text{mL}$)


The retention time is much smaller than I would have predicted from our calibration. 13.5 kDa should appear at ~200 mL based on that. Will need to look into this, but gel will confirm size tomorrow.

5/3/2017.

Gel of FL WT Pin1 is below. Definitely something fishy w/ the column. Peak is nearly identical.

UNICORN 5.20 (Build 500)
Result file: d:\...\dinusha\20170502 Pin1 WT Superdex 75 26 60001

- 20170502 Pin1 WT Superdex 75 26 60001:10 UV1 280nm
- 20170502 Pin1 WT Superdex 75 26 60001:10 UV2 260nm
- 20170502 Pin1 WT Superdex 75 26 60001:10 Cond
- 20170502 Pin1 WT Superdex 75 26 60001:10 Fractions



TITLE

PROJECT

DATE

Continued from page

I will fractions C1-C10. Adding 100 μ l of 10% NaN_3 .
A₂₈₀ is 0.198 @ 1mm. Conc is 94 μ M, w a volume of
 \sim 50 mL.

5/4/17 Mass Spec In Core Lab

Basic Protocol

Buffer A - 5% Acetic Acid in H_2O

Buffer B - Acetonitrile

Column: We've been using a Restek ¹⁵ 10 cm C18 column.


At NIH, they were using:

Temp = 40°C

- 1) Zorbax 300SB-C3 4.6x150 (PN 883995-909), or
- 2) Poroshell 300SB-C3 2.1x75 (PN 660750-909) (newer)

The C18 column probably doesn't resolve proteins as well as a C3 would, but our results have been okay.

Program



Continued to page

SIGNATURE

DATE

time (min)

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

DATE


Continued from page

Eric Anderson recommends ~ 0.2 mL/min, but we've had good results w/ ~~to~~ 0.5 mL/min. Pressure is ~ 100 bar (typically a bit less, goes down as Acetonitrile \uparrow)

Note that when the column is 1st ~~run~~ connected it should probably be run through a blank (inject buffer) or equilibrated for ~ 30 min. The blank can also help determine if anything is stuck on the column.

UV-Vis should be set to monitor 280 nm if proteins have Trp. 205 will catch proteins, too but it will catch buffers and small organics as well.

Vials are Sun-Sri (PN 200410). Make sure there are no air bubbles when loading vials. I typically load > 30 μ L into the vial and give them a good shake to eliminate air bubbles. These vials have a tapered insert and are disposable.



You may be able to ~~too~~ load less than 30 μ L, I haven't tried.

These settings are generally saved in the fitzkee protein method on the CPU desktop.

BUT you should always check everything to make sure.

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

DATE


Continued from page

LCMS Data for Pin1 CTD:


I injected 5 μ L of sample (approx ~133 μ M). This is density

(modified after loading)

VWD1 A, Wavelength=280 nm (D:\FITZKEE LAB\DATA\FITZKEE\INUSHA2017-05-0309-38-48\PROT000001.D)



MSD1 TIC, MS File (D:\FITZKEE LAB\DATA\FITZKEE\INUSHA2017-05-0309-38-48\PROT000001.D) ES-API, Pos, Scan, Frag: 70,



The column had contamination peaks near ~39 min. These diminished as we ran the column more, but and were present in the blank. The MS for these points was high MW but very messy, and I wonder if Colleen's polymers had contaminated the LC.

I intended to run 2 μ L, and I believe that would have been enough. I recommend 1 μ L when conc is ~100 μ M, as one can see contamination of CTD in the pin1 sample I ran next. (see p. 56)

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION

TITLE

PROJECT

DATE





Continued from page

The deconvolution is below. Expected mass is on page 43 and is 13157.0

Deconvolution Parameters

Adduct Ion(Positive): +H, 1.0079 Da
 Adduct Ion(Negative): -H, -1.0079 Da
 Low MW: 500
 DeconvStartChgMaximum Charge: 50
 Minimum Peaks in Set: 3
 Retain Residual: No
 Ion PWHH: 0.6 Da
 MW Agreement: 0.05 %
 Noise Cutoff: 1000 counts
 Abundance Cutoff: 10 %
 MW Assign: Curve fit
 MW Assign Cutoff: 40 %
 Envelope Cutoff: 50 %

Deconvolution of Spectrum # 1 @ 37.476 - 38.060 min



Component	Molecular Weight	Absolute Abundance	Relative Abundance
A	13155.73	133703	100.00
B	730.87	28281	21.15

TITLE

PROJECT


DATE

Continued from page


FL Pin1 construct: LCMS. Again, report says I injected 5 μ l, 1 would have been enough. Note the 2nd set of peaks at 39 min is much less.

(modified after loading)

VWD1 A, Wavelength=280 nm (D:\FITZKEE LAB\DATA\FITZKEE\INUSHA2017-05-0309-38-48\PROT000003.D)



MSD1 TIC, MS File (D:\FITZKEE LAB\DATA\FITZKEE\INUSHA2017-05-0309-38-48\PROT000003.D) ES-API, Pos, Scan, Frag: 70,



Note similarity between CTD and FL Pin1 (37.6 min vs. CTD 37.7)
This would be better on a C3 column I'm sure, as C3 would improve resolution.

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION




TITLE

PROJECT

DATE

Continued from page

Deconvolution of Spectrum # 1 @ 37.119 - 37.834 min



(FL Pin)

Component	Molecular Weight	Absolute Abundance	Relative Abundance
A	18578.15	604016	100.00

*** End of Report ***

Instrument 1 5/9/2017 6:38:45 PM Min

Page 1 of 1

Deconvolution parameters identical to what was done before. Expected MW is 18581.6 Both MW are slightly less than expected (2-3 Da). But the MS may need calibration, but I'm not concerned - the error is small enough that it isn't likely most likely not meaningful.

Continued to page

SIGNATURE

DATE

DISCLOSED TO AND UNDERSTOOD BY

DATE

PROPRIETARY INFORMATION