

Polarized Light




Optical Activity




From *Biophysical Chemistry, Part II* Cantor & Schimmel, Chapt. 8, p. 469

Basis of CD Spectroscopy



- Electric field dipole drives a current in the lock washer, which induces a magnetic dipole.
- Magnetic dipole must align with electric dipole $(\vec{\mu} \cdot \vec{m} \neq 0)$
- This is a significant oversimplification!

Basis of CD Spectroscopy



- If there's no mirror symmetry, E field doesn't induce a magnetic dipole.
- Main point: In chiral systems, electric fields can induce magnetic dipoles, which can interact with incident light itself

Secondary Structure

- Basis spectra based on model peptides (e.g. under helixforming conditions)
- Molar extinction coefficient per residue




Fig. 2. The CD for various secondary structures: α -helix (----), antiparallel β -sheet (-----), β -turn (----), and random coil (----), redrawn from Ref. 7.

DNA Structure




Figure 10.12 The CD of poly (dA) · poly (dT) as a function of temperature: 1.0°C (- - -); 38.8°C (-·-·); 44.7°C (--); 48.2°C (···); 58.3°C (--). [Adapted from data in J. Greve, M. F. Maestre, and A. Levin (1977), Biopolymers 16, 1489–1504. Reprinted from W. C. Johnson, "Circular Dichroism and Its Empirical Application to Biopolymers," in Methods of Biochemical Analysis, vol. 31, ed. D. Glick, copyright 1985 by John Wiley & Sons.]